

ADVANCED CARDIAC LIFE SUPPORT (ACLS) CART INVENTORY LIST

BACK OF CART - Backboard

SIDE OF CART

1. Handle Side
 - E-size O2 cylinder with regulator (replace tank if gauge reads < 1000 psi or < 1/2 full)
 - 1 - box gloves latex-free exam **medium (01194)**
2. Opposite Side
 - Suction yankauer tip w/ tubing **(00949)** (bottom rack) _____
 - portable suction unit with suction canister **(02797)**
 - replacement canister: suction canister 1200 ml **(00536)**
 - 1 - sharps container, medium **(Call CHS)**

TOP OF CART

1. 1" Notebook with:
 - Code Cart Check List - NIH Form #1355-1 (12-16, for ICU only) or 1355 (12-16)
 - ACLS Code Cart Inventory List (dated 12-16)
 - 3 ea - Emergency Response Records (NIH Form #2253)
 - Pediatric Emergency Drug Sheets – Generic (dated 6-16)
2. The following are **connected** and in good condition:
 - ECG cable
 - Quick Combo therapy cable
 - SPO₂ cable (with non-disposable sensor probe)
 - Defibrillation paddles present: (pediatric paddles for LIFEPAK 12 in bottom drawer)
 - NIBP cable (LIFEPAK 12 only)
3. **Supplies:**
 - 1 package - ECG radiolucent **electrodes (00516)** _____
 - 1 each: **Quick Combo** electrode: adult **(00498)**_____, pediatric **(00499)**_____
(2 ea. on ICU/Code Team carts only)
4. **Defibrillator:**
 - If LIFEPAK 20 Defibrillator:
 - 1 extra roll **ECG paper** - physio 50mm paper **(01632)**
 - Lifepak 12 Defibrillator (ICU/ Code Team/ PVCS/ OR/ PACU/ Cath Lab carts only)
 - 1 extra roll **ECG paper** - 100 mm/12 lead paper **(02943)**
 - Disposable **BP cuffs** - child **(01166)**, adult: small **(00492)**, average **(00491)**, large **(01165)**,
5. Code Team Cart Only –
 - Keys: Elevator override key, AED box key
 - iSTAT, portable **ETCo2** (maintained by Respiratory Therapy Section)
 - 12 lead-compatible ECG Cable (instead of 5-lead cable)
 - Disposable **BP cuff** - infant **(01167)** - Code Team cart only)
 - GlideScope Ranger and associated ET tubes for adults and pediatric patients

FRONT OF CART - DRAWER 1: Medication Tray (contents sealed in plastic bag).

FRONT OF CART

DRAWER 2: Needles, Syringes

<p>2- spinal needles: 20ga. x 3½" (01121) ____ 1- 60 ml Leur (01073) ____ 1- 30 ml Leur (01072) ____</p>	<p>Protected IV catheters:</p> <p style="padding-left: 40px;">2 - 16 ga.(01142) 2 - 18 ga.(01143) 2 - 20 ga.(01144) 2 - 22 ga.(01146) 2 - 24 ga (01147)</p> <p>IV Start Kit: <i>*Place one of each item in clear ziploc bag (00633) and label "IV Start"</i></p> <ul style="list-style-type: none"> • 2 - 2x2 gauze (00715) ____ • 2 - DSG IV 3000, 4"x5" (00987) ____ • 2 - Luer 7" extension set (03281) • 2 - tourniquets - latex free (01964) • 2 - Chlorhexidine 1.5 ml swabs (02708) ____ • 2 - VenoGuard (00773) ____ • 1 - NS 10 ml pre-filled syringes (order from pharmacy) ____ 		<p>Messenger Kit <i>*Place contents of messenger kit in ZipLoc bag (00633)</i></p> <ul style="list-style-type: none"> • 1 each: pen, magic marker • 1 - 6 x 9 lab STAT specimen bags (02326) <ul style="list-style-type: none"> * 1 - Clinical Laboratory Ordering Record CHEMISTRY (NIH 2901 (4-11) * HEMATOLOGY (NIH-2902 (2-11) • 3 - 3x5 cards w/ CODE BLUE label affixed • 3 - extra CODE BLUE labels • 3 - plain white labels • 6 - MEDICATION ADDED labels (03052) <p>4 ea: gauze 2x2 dressing (00715) ____</p> <p>2 - Chlorhexidine swabstick (02599) ____</p>
<p>5 - Normal Saline 10 ml pre-filled syringes (obtain from pharmacy) ____</p>			
<p>8 - 10ml luer w/o needle (01068)</p>			
<p>3 - 3 ml luer w/o needle (01060) 2 - 3ml with needle (02447)</p>			
<p>3 - 1ml (01890) 2 - 1ml with needle (02456)</p>	<p>2 - CTSD microclave male plug (03186)</p> <p>2 each - Butterfly needle, safety:</p> <ul style="list-style-type: none"> • 21 ga.(02488) • 23 ga.(02487) 	<p>3 - needle Safety Glide: ____ 25 ga. x 5/8" (02454)</p> <p>5- needle (unprotected) 19 ga x 1 1/2" (00825)</p> <p>3 - blunt fill needle-filter: 18 ga. (02754) ____</p>	<p>10 - vial access cannula (01094) 2 - BD Q-Syte Multi-Dose Adapter (03279)</p> <p>1 - 2" bandage, EBAN colored (00186) 1 - 1" water proof cloth tape (00736) 1 - 1" roll transpore tape (01337) 5 - adhesive bandages (band-aids) (00016) 15 - alcohol swab pads (00718)</p>
<p>2 - stopcock, double 4-way (00049) ____ 2 - stopcock, 4-way lipid resist. (01043) ____</p>			

DRAWER 3: IV fluids / Tubing

- 2 - sol NS: 1000 ml **(00701)** _____
- 1 - sol NS: 500 ml **(00702)** _____
- 1 - sol D₅W: 250 ml **(00939)** _____
- 1 - sol NS: 250 ml **(00940)** _____
- 2 - sol D5W: 100 ml **(00696)** _____

1 - Dopamine PREMIX (400mg/250 ml D₅W) **(Order from: Hospital Services - Visual Supply Catalog/ Pharmacy Emergency Stock/ Code Cart Replacement Meds)**

- 2 - set IV admin. (primary) Carefusion **(03278)** _____
- 1 - set IV admin. universal(secondary) Alaris **(00493)** _____
- 1 - IV tubing, low-sorbing Carefusion set **(03263)** _____

- 1 - set extension tubing 20" **(03151)**
- 1 - armboard IV supporter **(00613)**
- 1 - infusion pressure bag **(01042)**

DRAWER 4: Dressings/Tubes

- 2 - underpad (chux) **(00109)**
- 4 - mask, Technol procedure **(00244)**
- 2 - glasses, protective eyewear **(01307)**

- Gloves, surgical: 1 each size
- 6 **(02811)** _____, 6 ½ **(02812)** _____
 - 7 **(02813)** _____, 7 ½ **(02814)** _____
 - 8 **(02815)** _____, 8 ½ **(02816)** _____

1 - feeding tube 8 Fr. **(01216)**

Salem Sump tubes:

- 1 - 12 fr. **(02371)**
- 1 - 16 fr. **(00778)**
- 1 - syringe, catheter tip 60 ml **(01074)** _____
- 3 - Surgilube packets **(from Pharmacy)**
- 1 - hand sanitizer, pocket size **(01657)** _____
- 1 - NS 0.9% Irrigation - 500 ml bottle **(00960)** _____
- 2 - dressing (combine 8x7) **(00734)** _____

- 6 - gauze 4x4 dressing **(00743)** _____
- 2 - disposable cold packs **(02132)**

- 1 - scalpel w/handles, sterile: #11 **(00784)** _____

- 1 - scissors, sterile 5½" blunt **(00789)** _____
- 1 - Kelly forceps straight, sterile **(01217)** _____
- 2 - suture: 2-0 silk black w/cutting needle **(01040)**

- 2 - Benzoin swab sticks **(00020)** _____
- 1 - Povidone iodine, 2 oz. bottle **(00157)** _____

CT/ICU Carts only - Substitute

- 1 - Povidone iodine, 8 oz. bottle **(01654)** _____

DRAWER 5: Respiratory

O2 Delivery/ Suction

1 ea - Non-Rebreather Mask: adult & pediatric & infant **(infant on ICU CARTS ONLY)**

1 ea - Nasal Cannula:- adult, pediatric, & infant

O₂ Flowmeter w/ green adapter

1 ea - Nebulizer Kit

(Contact RT for above)

ABG Kits: 2 ea. **(02586)** _____

Yankauer suction device - 1 ea:
adult **(00950)**, pediatric **(01582)**

Suction catheter kits - 2 each (Fr.): 10 **(00233)**,
12 **(02538)**, 14 **(00167)**

1 each ICU CARTS ONLY:

6 Fr (02156); 8 Fr (01700)

1 ea - Aerosol Face Mask (Code Team cart only)

Bag/ Mask

Resuscitation Bag/Valve/Mask:

1 ea: adult & pediatric & infant **(infant on ICU CARTS ONLY)**

Mask, soft seal clear: 1 ea: Adult size 6, adult size 5, toddler, child, infant, & neonate **(neonate on ICU CARTS ONLY)**

(Contact RT for ALL MASKS above, check expiration date)

Laerdal Pocket Mask: 1 ea **(00204)**

Carts with LIFEPAK 12 - add:

1 ea - adult disposable pulse oximeter (OxiMax) sensor **(01227)**

1 ea - pediatric disposable pulse oximeter (OxiMax-P) sensor **(01041)**

** Code Team only: 1 ea.- infant disposable pulse oximeter (OxiMax-I) sensor **(00286)***

Intubation Roll

MACGILL FORCEPS - 1 ea: Large, small **(Call CHS)**

STYLETS - 2 ea: 14 Fr. **(00078)**, 6 Fr. **(call CHS)**

NASOPHARYNGEAL AIRWAYS - 1 ea:

Adult: 6 **(01606)**, 7 **(01554)**, 8 **(01553)**, 9 **(02621)**

Pediatric (red rubber): 16 Fr **(01124)**, 18 Fr **(01125)**

Benzoïn swabstick: 1 ea **(00020)** _____

Tongue blades: 1 ea. **(00033)**

1" water proof cloth tape roll: 1 ea. **(00736)**

Syringe, 10 ml slip tip: 1 ea. **(01067)**

Laryngoscope HANDLES: 1ea – medium/ small

*Batteries for Handles – 2 each:

Size C **(01501)** _____

Size AA (02936) _____

Laryngoscope BLADES: 1 ea- Mac (curved): #1,#2,#3,#4

Miller (straight): #0, #1, #2, #3 **(Call CHS)**

*spare laryngoscope bulb: 1 ea. large **(02331)**, small **(02285)**

*Code Team / ICU ONLY: Disposable Laryngoscope **(Call CHS)**
Greenline Laryngoscope HANDLES (non-disposable) - Adult,*

Pediatric

ENDOTRACHEAL TUBES - 1 each:

Uncuffed- 3 **(01706)**____, 3.5**(01778)**____4**(01779)**____,
4.5**(01780)**____ 5 **(01781)**____, 5.5**(01782)**____
6 **(01783)**____

Cuffed – 4 **(01786)**____, 4.5**(01787)**____ 5**(01788)**____,
5.5**(01789)**____ 6 **(02384)**____, 6.5**(02385)**____
7 **(02386)**____, 7.5**(02387)**____ 8 **(02388)**____
8.5**(02389)**____ 9 **(02390)**____

OROPHARYNGEAL AIRWAYS (1 ea): 5 **(02637)**, 6 **(02638)**,
7 **(02639)**, 8 **(00201)**, 9**(00202)**, 10 **(01696)**

Lidocaine Jelly tube: 1 ea. **(from Pharmacy)** _____

Detectors

ETCO2 detectors (Easy-cap) - 1 each:

		adult (02357)____, pediatric (02356)____ *ICU and Code Team Cart only: 2 each: CO2 detectors
--	--	---

Drawer 6 (BOTTOM): Lines/

<p>Defibrillation gel: 1 tube (00531)</p> <p>Flashlight: 1 ea (00681) with 2 extra - "D" batteries (02312) _____</p> <p>Stethoscope, dual head: 1 ea. (00578)</p> <p>BP cuffs multi-complete set (small, adult, large, child) (02840) with sphygmomanometer attached.</p> <p><u>Code Team Cart Only:</u> add 1 ea: Infant BP cuff (special order - CHS) Heavy duty trash bag Towel Roll</p>	<p>Lines</p> <p>Triple Lumen Catheter Kits: 1 ea. 7 Fr. x 16 cm (00285)____ 5.5 Fr. x 8 cm Pedi (00279) _____ <i>(CT & ICU Carts Only – 2 ea.)</i></p>	<p>Code Team & ICU Carts Only:</p> <p>2 - Sternal - Iliac Needle (Intraosseous) with adjustable length 3/16" – 1/8" (00243) _____</p> <p>1 - pneumothorax kit (00297) _____</p> <p>1 - cricothyrotomy kit (00302) _____</p> <p>2 each: non-protected angiocath: 14 ga. x 2"(02966) _____ 16 ga. x 2"(02967) _____</p> <p>Disposable laryngeal mask airways (LMA): 1 each #1 (02915) _____, #1.5 (02916) _____ #2 (02917) _____, #2.5 (02918) _____ #3 (02919) _____, #4 (02920) _____ #5 (02921) _____</p> <p>Defibrillation gel pads – 4"x9"(01674) _____ 6"x9"(01675) _____</p> <p>Pocket calculator Pen</p> <p>ICU CART ONLY: EZ-IO drill and needles (one of each: 15, 25, 45mm). ASK RESPIRATORY TO REPLACE ITEMS RELATED TO THE EZ-IO. PLEASE REMOVE THIS ITEM FROM THE CART BEFORE SENDING IT BACK TO CHS FOR RESTOCKING.</p>
--	---	--

**** IMPORTANT**** Prior to locking the cart assemble the following:

- stethoscope

- Adult BP cuff
- manual resuscitator bag-mask device

Check for proper functioning of laryngoscope handles and blades (batteries strong and lights shine brightly on each blade).

****Highlighted items expire – check expiration dates**